

Ordering Information for Arrhythmogenic Right Ventricular Dysplasia/Cardiomyopathy (ARVD/C) Testing

Indications for Testing

- Clinical diagnosis of Arrhythmogenic Right Ventricular Dysplasia/ Cardiomyopathy (ARVD/C)
- Unexplained cardiovascular symptoms such as palpitations and/or history of syncopal episodes in young adults, especially athletes
- Family history of ARVD/C
- Family history of sudden cardiac death in individuals under age 45

Ordering Information for Single Gene Tests

Gene(s)	CPT Codes	Test Code
PKP2	83891(1) 83892(1) 83898(15) 83904(30) 83909(30) 83912(1)	190701
DSP	83891(1) 83892(1) 83898(38) 83904(76) 83909(76) 83912(1)	190702
DSC2	83891(1) 83892(1) 83898(17) 83904(34) 83909(34) 83912(1)	190703
DSG2	83891(1) 83892(1) 83898(17) 83904(34) 83909(34) 83912(1)	190704
TMEM43	83891(1) 83892(1) 83898(12) 83904(24) 83909(24) 83912(1)	190705

Ordering Information for Multi-Gene Panels*

PKP2, DSP, DSC2, DSG2, TMEM43	83891(1) 83892(1) 83898(99) 83904(198) 83909(198) 83912(5)	190799
--------------------------------------	---	--------

Family Testing (single amplicon)

(applies to all genes)	83891(1) 83892(1) 83898(1) 83904(2) 83909(2) 83912(1)	use single-gene test code
-------------------------------	--	---------------------------

Test Methodology and Sample Requirements

- Amplification by polymerase chain reaction (PCR); sequencing of entire protein-coding region
- For blood samples:
 - 2 mL whole blood in EDTA tube (lavender top)
 - Samples can be stored briefly at 4°C, but should be shipped on day of collection.
- For buccal swab samples: (only accepted for family testing)
 - Please contact client services at 1-866-647-0735 for instructions.
- All sample types should be shipped overnight at room temperature.
- To request a sample shipping kit, please call 1-866-647-0735.

For more information, please contact Correlagen Diagnostics, Inc., at 1-866-647-0735 or visit us on the web at www.correlagen.com.